# “FIN ARTHROPROSTHESIS”

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**Disegno distale**  
*Distal design*

Advanced biomechanical concepts, cutting-edge design, careful selection of materials and high quality technology have been harnessed to meet the needs of modern orthopaedics. And this brings undisputed benefits for patients.

**Testa protesica  
in Biolox®Delta  
Biolox®Delta  
Prosthetic head**

## TECHNICAL AND GEOMETRICAL DESCRIPTION

The "Fin" stem is a straight stem designed for easy insertion and efficacious anchorage.

### Features

#### A) PROXIMAL DESIGN

- Dorsal fin ensuring rotational stability.
- Wedge shape, with the longest side towards the trochanter, allowing press-fitting of the stem.

- Small support collar on the calcar to ensure axial stability and prevent re-absorption.
- Longitudinal normalization lines specially designed for physiological distribution of all forces.

#### B) CENTRAL DESIGN

- The middle part is appropriately sized to allow a gradual reduction of stresses transferred to the bone from the proximal to the distal zone.

#### C) DISTAL DESIGN

- The tapered distal part prevents any contact between the stem and the cortical walls thus avoiding distal interposition of fibrous matter.

## MATERIALS:

### BASE MATERIALS

The prosthesis is made of titanium with a low elasticity modulus (can tolerate temporary deformation under stress) to encourage balanced load distribution over all the surrounding areas. No concentrated forces are applied to any particular area or element, as in a healthy, unbroken femur.

The titanium alloy used has the following characteristics:

- high purity achieved by means of the ELI (Extra Low Interstitial) process to increase resistance to the propagation of fatigue cracks (ELI grade 5)
- optimum biocompatibility
- optimum resistance to corrosion
- optimum mechanical strength (forged)
- conforms to American (ASTM) and British (BS) international standards

### COATING MATERIALS:

The Y 370 VPS porous titanium coating permits mechanical anchorage between the prosthesis and the bone without the interposition of fibrous tissue.

To encourage bone ingrowth the proximal part of the FIN stem is coated with oxide-free, 99.4% pure titanium (Eurocoating® Y 370 VPS) using a vacuum plasma spray technique.

The main characteristics of the coating are the following:

- No variation in the mechanical strength of the base material thanks to low temperature application of the coating ( $\pm 300^{\circ}\text{C}$ ).
- Excellent anchorage (greater than 40 MPA) between the coating and the titanium alloy.
- Porosity in the range 20-65% with diameters of 100-250 microns.
- High stability coating particles.
- 300 micron thickness.

### HEAD MATERIALS:

- Bioceramic BIOLOX®FORTE, BIOLOX®DELTA
- Chrome and Cobalt Alloy
- Stainless steel with high nitrogen content

## NOTES:

- The "Fin" stem is compatible with the biomechanics of the articulation: ccd angle  $135^{\circ}$ .
- In line with the European standard 12-14 taper ( $5^{\circ} 42' 30''$ ).

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
Fig. 1

## ADVANCED TECHNOLOGY COATING:

- **Proximal part.** The porous Y 370 VPS coating in the proximal area encourages bone ingrowth. The photographs alongside illustrate two metallographic sections of the Y 370 VPS coating, 300 microns thick (see figures 1 and 2).
- **Central part.** The central part is characterized by its corundum surface treatment which contributes to the stem's secondary stability.
- **Distal part.** The reduced porosity of the distal part (achieved by glass bead surface treatment) allows most of the stresses to be transferred proximally.


Fig. 2

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


## STABILITY:

### a) PRIMARY:

- **Dorsal proximal fin**  
Guarantees rotational stability. (Detail A)
- **Small proximal collar**  
Guarantees axial stability and prevents longitudinal subsidence. (Detail B)
- **Proximal press-fit geometry**  
Guarantees excellent proximal metaphyseal filling. (Detail C)

### b) SECONDARY:

Secondary stability is ensured by coating the proximal part in Y 370 VPS. This encourages bone ingrowth and can be seen in details A and B.

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## SIMPLIFIED IMPLANT TECHNIQUE:

Thanks to the shape of the stem and the use of a minimum number of universal instruments.

## SIZING:

To achieve optimum results with high numbers of patients the "Fin" stem is available in seven different lengths, varying proportionally to diameter, to cover a wide range of femur sizes. Available sizes: 8 - 9 - 10 - 11 - 12 - 13 - 14 - 15 and 16 mm

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
### TESTINE PROTESICHE IN CERAMICA - CERAMIC PROSTHETIC HEADS

Codice Code	Cono Taper	Collo Neck	C.I.R. * R.I.C. *
11-0230105.E	12/14	corto/short	-3,5 mm
11-0230110.E	12/14	medio/medium	0,0 mm
11-0230115.E	12/14	lungo/long	+3,5 mm

ø **28** mm e (materiale Al<sub>2</sub>O<sub>3</sub>) BIOLOX®FORTE  
ø mm **28** and (Al<sub>2</sub>O<sub>3</sub> material) BIOLOX®FORTE

Codice Code	Cono Taper	Collo Neck	C.I.R. * R.I.C. *
11-0240105.E	12/14	corto/short	-3,5 mm
11-0240110.E	12/14	medio/medium	0,0 mm
11-0240115.E	12/14	lungo/long	+3,5 mm

ø **32** mm e (materiale Al<sub>2</sub>O<sub>3</sub>) BIOLOX®FORTE  
ø mm **32** and (Al<sub>2</sub>O<sub>3</sub> material) BIOLOX®FORTE

Codice Code	Cono Taper	Collo Neck	C.I.R. * R.I.C. *
11-0240135	12/14	corto/short	-4,0 mm
11-0240140	12/14	medio/medium	0,0 mm
11-0240145	12/14	lungo/long	+4,0 mm

ø **36** mm e (materiale Al<sub>2</sub>O<sub>3</sub>) BIOLOX®FORTE  
ø mm **36** and (Al<sub>2</sub>O<sub>3</sub> material) BIOLOX®FORTE

TESTINE PROTESICHE IN CERAMICA - CERAMIC PROSTHETIC HEADS

Codice Code	Cono Taper	Collo Neck	C.I.R. * R.I.C. *
11-0240205	12/14	corto/ <i>short</i>	-3,5 mm
11-0240210	12/14	medio/ <i>medium</i>	0,0 mm
11-0240215	12/14	lungo/ <i>long</i>	+3,5 mm

ø 28 mm BIOLOX®DELTA  
ø mm 28 BIOLOX®DELTA

Codice Code	Cono Taper	Collo Neck	C.I.R. * R.I.C. *
11-0240305	12/14	corto/ <i>short</i>	-3,5 mm
11-0240310	12/14	medio/ <i>medium</i>	0,0 mm
11-0240315	12/14	lungo/ <i>long</i>	+3,5 mm
11-0240320	12/14	extra lungo/ <i>extra long</i>	+7,0 mm

ø 32 BIOLOX®DELTA  
ø mm 32 BIOLOX®DELTA

Codice Code	Cono Taper	Collo Neck	C.I.R. * R.I.C. *
11-0240405	12/14	corto/ <i>short</i>	-4,0 mm
11-0240410	12/14	medio/ <i>medium</i>	0,0 mm
11-0240415	12/14	lungo/ <i>long</i>	+4,0 mm
11-0240420	12/14	extra lungo/ <i>extra long</i>	+8,0 mm

ø 36 mm BIOLOX®DELTA  
ø mm 36 BIOLOX®DELTA

TESTINE PROTESICHE IN ACCIAIO - PROSTHETIC HEADS IN STAINLESS STEEL

Codice Code	Cono Taper	Collo Neck	C.I.R. * R.I.C. *
11-0205105.E	12/14	corto/ <i>short</i>	-3,5 mm
11-0205110.E	12/14	medio/ <i>medium</i>	0,0 mm
11-0205115.E	12/14	lungo/ <i>long</i>	+3,5 mm
11-0205120.E	12/14	extra lungo/ <i>extra long</i>	+7,0 mm
11-0205125.E	12/14	xx-lungo/ <i>xx-long</i>	+10,0 mm

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The cemented "Fin" stem is a straight stem with a collar and a geometrical shape that simplifies the implant technique.

### Features

#### A) PROXIMAL DESIGN

- **Dorsal fin** ensuring rotational stability and reducing stresses on the cement.
- **Wedge shape**, allowing greater metaphyseal filling and reducing the cement "coat" to a minimum.
- **Small support collar** on the calcar to ensure axial stability and compression of the cement.
- **Longitudinal normalization lines** specially designed for physiological distribution of all forces.

#### B) CENTRAL DESIGN

- The middle part is appropriately sized to allow a gradual reduction of stresses transferred to the bone-cement interface from the proximal zone.

#### C) DISTAL DESIGN

- The tapered distal part prevents any contact between the stem and the bone.

## MATERIALS:

The material used for the cemented "Fin" stem is high nitrogen content stainless steel, to ISO 5832/9, the best material, in the present state of the art, for cemented prostheses. Thanks to the high nitrogen content and the final post-casting finish, it represents a marked improvement on the characteristics of the stainless steel most widely used in implant surgery - ISO 5832/1. ISO 5832/9 stainless steel is characterized by:

- High resistance to corrosion.
- Excellent mechanical strength.
- Responds better to cold hardening, for high resistance to fatigue.

## SIZING:

To achieve optimum results with high numbers of patients the "Fin" stem is available in nine different lengths, varying proportionally to diameter, to cover a wide range of femur sizes and reduce the amount of cement needed to a minimum.

Available sizes: 8 - 9 - 10 - 11 - 12 - 13 - 14 - 15 - 16 mm

Lengths: 142 - 146 - 152 - 158 - 162 - 167 - 170 - 172 (L 220) - 177 (L 220) mm.

## EXTERNAL FINISH:

The external part is bombarded with ceramic microspheres which increase the already high mechanical surface strength and create a micro roughness that guarantees a perfect stem-cement interface.

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The cemented Fin stem can be implanted using the same instruments as for the cementless stem but rasping one size bigger than the one desired.

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**11** **Manico a "T" per alesatore conico**  
*"T" handle for conical reamers*

Cod. 11-0381520

**12** **Attacco rapido jacobson**  
*Quick coupling for conical reamers*

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GRUPPO BIOIMPIANTI S.r.l.